

PRINT DATE: Tuesday, 12 March 2024



► VIEW DATASHEET ON-LINE [www.chinesport.com/go/AR10044/](http://www.chinesport.com/go/AR10044/)

#### AR10044 SMALL BASE WITH HOLES

It's a base module of small rectangular shape, with thirty holes to insert some specific elements freely. This base module used to attach different aids is made of highly resistant, anti-bacterial material of the highest quality, to last over time without alterations; in particular, this base is resistant to water and UV rays and can be easily disinfected.

Several work tools have been conceived to make exercises of grasping and inserting in the holes of this base module: 1) cylinders and, 2) spheres, 3) pins with magnet and colored discs.

The cylinders are of two types differing for diameter and material, instead the spheres can be of two different diameter. Lastly the pins with magnet, when inserted in the base, offer the additional exercise of applying and removing special coloured discs supplied as standard with these pins. Moreover as alternative to the coloured discs just mentioned it's possible to use discs with white surface to write on numbers and letters for instance.

The exercises range can be extended with grasping and inserting additional cylinders one to another as well as special big disc can be inserted in the cylinders or used for paths passing through the inserted cylinders.

Another variety of exercises can follow the application of cylinders on the module base such as grasping and inserting a ring attached to an elastic. So the elastic can also make a track passing through the cylinders in before ending the exercise with the ring inserted in the final cylinder.

This module can be easily carried by the therapist or assistant in particular for home therapy. In some countries home therapy is well recognized and sponsored to avoiding hospital queues and for a more patient involvement and constant commitment.

Dimensions: 39.5 cm x 17.5 x 1.7 h cm; Weight: 0.9 kg



AR10044 - HOLE BASE - SMALL MODEL

## Technical specifications and dimensions

Dimensions : 39.5 x 17.5 x 1.7 h cm

Set weight : 0.9 kg

## Documentation available online

<b>DICO</b> - Occupational therapy English		.PDF 1.5 MB
<b>DICO</b> - Terapia occupazionale Italian		.PDF 1.5 MB
<b>IFU</b> - Configurable aids English		.PDF 4.4 MB
<b>IFU</b> - Sussidi configurabili Italian		.PDF 4.4 MB
<b>Manual</b> - User and maintenance manual Italian, English		.PDF 249.24 KB

 The document is ready for download in reserved area.

**Accessories**



**AR10139 FANTASY STICKS SET 5**

A series of 8 aluminium cylinders with 15 mm diameter. This type of cylinder offers the option of inserting them one on top of the other.



**AR10141 FANTASY STICKS SET 6**

A series of 8 plastic cylinders. This type of cylinder offers the option of inserting them one on top of the other.



**AR10034 OLYMPIC DISCS SET 2**

A series of 8 discs in four colours. The half-moon discs have a grip that is easier to grasp, using the hand in various positions.



**AR10143 SPHERES SET Ø20/Ø40**

A series of 8 spheres made from plastic material: 4 spheres with Ø 20 mm and 4 with Ø 40 mm that can be inserted into a prearranged base module.



**AR10135 MAGNETIC GAMES SET**

A series of 8 pins with a magnet at the top and 8 metal discs in 4 colours.



**AR10036 MAGNETIC GAMES SET 2**

A series of 16 half-moon metal discs in 4 colours. The grasp can be terminal digital or digito-lateral, involving the distal phalangeal and wrist joints.



**AR10037 MAGNETIC GAMES SET 3**

Series of ten metal discs with a special, rewritable white film. They can be used to write numbers and letters individually, creating a memory game.



**AR10123 ELASTIC SLALOM**

Made to exercise the subject in "slalom" tasks between the pegs. It requires precision in movement, the ring facilitates grip.

**More pictures of product**



AR10044 - HOLE BASE - SMALL MODEL

