

PRINT DATE: Tuesday, 12 March 2024



► VIEW DATASHEET ON-LINE www.chinesport.com/go/AR10048/

AR10048 TRACKS - LARGE BOARD

This aid comes with two different sized aluminium pins, two different sized spheres, and one plastic cylinder as in item AR10047, however the base module has a larger square shape, with several patterns or possible paths.

The elements offer various grips and types of "grasp".

It improves the use of fingers and wrists, both in flexo-extension and pronosupination, and improves the use of elbows in flexo-extension, hand-eye coordination both in adult and adolescent patients who need to learn and practise skills they have never used before or have lost after a trauma injury to the limbs or due to neurological lesions.

Helps improve writing skills and using keys in locks.



GRASP AND INSERT



USE WITH ONE HAND



USE WITH TWO HANDS



HORIZONTAL AID POSITIONING



VERTICAL AID POSITIONING



AR10048 - TRACKS - LARGE MODEL

Technical specifications and dimensions

Dimensions : 40.5 x 38 x 6 h cm

Weight : 2.4 kg

Documentation available online

DICO - Terapia occupazionale Italian		.PDF 1.5 MB
DICO - Occupational therapy English		.PDF 1.5 MB
IFU - Sussidi per inserire ed afferrare Italian		.PDF 3.28 MB
IFU - Grip and insert exercises English		.PDF 3.27 MB
Compliance statement - Declaration of conformity Italian, English, German, French, Spanish		.PDF 123.46 KB
Manual - User and maintenance manual Italian, English		.PDF 256.32 KB

 The document is ready for download in reserved area.

More pictures of product



AR10048 - TRACKS - LARGE MODEL

